

CENTRACARE Laboratory Services

MEMO

TO: CCLS Clients
FROM: CCLS - Specimen Referral Center
DATE: 6/10/2019
SUBJECT: Histamine, Plasma

Effective immediately, Mayo test ID **FHSPL** will reflect the following changes.

HISTAMINE, PLASMA (LAB0234369)

Explanation: The specimen stability and analytic time is changed.

| Current Specimen Stability | | |
|----------------------------|-------------|----------|
| Specimen | Temperature | Time |
| Plasma EDTA | Frozen | 365 days |

| New Specimen Stability | | |
|------------------------|-------------|---------|
| Specimen | Temperature | Time |
| Plasma EDTA | Frozen | 28 days |

| Current Analytic Time |
|-----------------------|
| 2 days |

| New Analytic Time |
|-------------------|
| 2-3 days |

Please feel free to contact the CCLS Specimen Referral Center with any questions at 320-251-2700 Ext. 57320